




**BARNSTEAD
B-PURE**



**Thermo Scientific
Barnstead Water Purification
Cartridge and Filter Systems**



exceptional value

versatile configurations

Thermo
SCIENTIFIC

Thermo Scientific Barnstead Cartridge and Filter Systems

Achieve the desired water types for your specific applications

Our complete line of cartridge and filter systems offers solutions for single step and pretreatment purification at an exceptional value. Cartridge systems are excellent for use with humidification/environmental chambers, autoclaves, water baths, CO₂ incubators, stills, ultrapure water systems, and eliminating corrosion and scaling when used with cooling loops for lasers and radar systems.



Versatile Configurations

- Versatile configurations for an array of specific applications
- Customize modular systems to remove specific impurities. Over-purifying water can waste time and resources, not to mention the possibility of being detrimental to your application
- Quick-release designs provide easy cartridge change-outs and simplify maintenance, bypassing the costs of service calls
- Wall-mounted systems free up bench space
- Cartridge systems are available for installation in pressurized and non-pressurized water line environments
- Wide selection of sizes supports space limitations and volume needs






Exceptional Value

- Robust polypropylene cartridge casings lengthen system life and prevent impurities leaching into the purified water
- High quality resins reliably remove targeted impurities
- Purity (resistivity) meter options include Pura-Lite “Go/No Go” indicators and digital purity meters for quick-glance purity values and signs cartridge replacement is needed
- Instead of joints and seams that are chemically adhered, ours are ultrasonically welded to prevent contamination and allow for recycling



Thermo Scientific Barnstead Cartridge and Filter Systems Quick Select Guide



Pressurized Cartridge Systems



	B-Pure™ Single Holder	B-Pure Double Holder	B-Pure Half-size Holder	High Temperature Stainless Steel	Economy Filter Holder	Harvey DI+
Digital Purity Meter		💧				
Analog Purity Meter	💧	💧		💧		💧
Color Change Indicator						
Remote Dispenser						💧
Wall-Mounted	💧	💧	💧	💧		💧
High Temperature Tolerant				💧		
Quick-Release Canisters	💧	💧	💧			💧
Maximum Flow Rate L/min	4	4	2	3	2	4

* All Hose Nipple Cartridges except 7373D contain dyes which change color when resin is exhausted

Identify the best system based on your unique requirements

- Will the cartridge system be installed in a pressurized or non-pressurized water line?
- What is your principal application for pure water?
- What do you need removed from your water?
- What is the quality of the feed water?
- What is the water volume you need on a daily basis?
- Do you have temperature or flow rate parameters?

Non-pressurized Cartridge Systems



Hose Nipple*	Bantam™ Deionizer	Benefits	Typical Uses
		See real-time purity values	With more sensitive applications
	💧	Indicates when to install new cartridge at quick glance	If cartridge system feeds equipment, such as autoclave, where real-time resistivity is not required
💧		Indicates when to install new cartridge	When feeding lab equipment
		Hand-dispense water from the cartridge system	For rinsing glassware or filling water bath
💧	💧	Saves laboratory bench space	For autoclaves and dishwashers
		For use where water temperature exceeds 37.8°C (100°F)	Good option for cooling loops
		Fast and simple cartridge change-outs	For all applications
1.25	0.63		

What You Need to Remove from Feed Water Determines the Right Cartridge

Organics or Chlorine


Organic solids are present in water sources from plant and animal decay as well as human activity. These may include proteins, alcohols, chloramines, and residues of pesticides, herbicides, and detergents which foul ion exchange resins and interfere with organic analysis. Chlorine is an additive used to treat drinking water and may need to be removed for biological protocols. Activated carbon will remove both organics and chlorine from the feed water.

Resin needed - activated carbon

Colloids

Colloidal particles typically have a slightly net negative charge ranging in size from 0.01 to 1.0µm and can be either organic or inorganic. Unlike suspended particles, colloids do not settle out by gravity, but remain suspended in liquid. Colloids clog filters, interfere with instrument operation, foul reverse osmosis membranes, and can bypass ion exchange resins, resulting in lower resistivity in deionized water systems. Macroreticular resins enhance performance of deionization resins and extend final filter life.

Resin needed - macroreticular resin



Feed water with organic and chlorine contaminants come into contact with the activated carbon in the cartridge. The impurities adhere to the surface of the carbon, purifying water.



Inorganic Ions

Impurities such as silicates, chlorides, fluorides, bicarbonates, sulfates, phosphates, nitrates, and ferrous compounds are present as cations (positively charged ions) and anions (negatively charged ions). Water with a high concentration of ions will conduct electricity and have high conductivity and low resistivity, as the two are inversely related. Ions will adversely affect the results of inorganic analyses such as IC, AA, ICP/MS, and may retard cell and tissue growth in biological research. These ions can also affect the cartridge life in deionized water systems.

Resin needed - ion exchange resin, cation exchange, or anion exchange

- Two-Bed (high capacity) resin – contains a layer of cation resin followed by a layer of anion resin within the same cartridge. This resin will not remove ions as thoroughly as the mixed-bed, but it will last longer.
- Mixed-Bed (ultrapure) resin – contains a mixture of anion and cation resins to improve ion exchange reaction kinetics. It produces high resistivity water at a pH of 6.8 to 7.2 that is free of silica, carbon dioxide, and other ionized impurities.

Oxygen

To help prevent corrosion in applications such as closed cooling loop systems for radar, computer, and laser systems choose an oxygen removal resin.

Resin needed - oxygen removal



Only high quality resins are used in our cartridge and filter systems to ensure consistently high quality product water.

Thermo Scientific Barnstead Pressurized cartridge systems



Digital Purity Meter provides quick-glance notification of when to install a new cartridge

B-Pure Cartridge Holder Systems

These expandable, modular systems are the best choice when applications require up to 30 L/day with a maximum flow rate of 4 L/min. B-Pure systems are engineered to support feed water pressures up to 100 psi and feature a simple quarter turn, quick-release cartridge exchange design. Standard wall-mounting brackets are included. Interconnect B-Pure Holders to form custom water treatment systems.

- **B-Pure Single Holder** produces water at up to 4 L/min and can be used separately or in series with other single holders. Pura-Lite "Go/No Go" analog resistivity indicators are available.
- **B-Pure Double Holder** includes 2 canisters and draw-off faucet. Water is dispensed at 4 L/min and optional digital or analog purity meters are available.
- **B-Pure Half-size Holder** supports flow rates up to 2 L/min with an optional pressure gauge kit to monitor inlet and outlet water pressure. Can support either half-size cartridges or the 10 in (25.43 cm) filters.

Typical applications - pretreatment, deionization, filtration

B-Pure systems offer easy-to-use, easy-to-maintain cartridge replacement

Step 1:
Change cartridge when the resistivity is below the desired level or the cartridge is more than 1 year old. For filters, you'll notice a significant decrease in pressure or flow rate

Step 2:
Remove the canister from the head by depressing the thumb lever and rotating the hand ring quarter turn to the left

Step 3:
Remove exhausted cartridge and discard

Step 4:
Replace old cartridge with a new one

Step 5:
Install canister by depressing thumb lever and rotating quarter turn to the right, until the locking pin is in the appropriate position

Other Pressurized Cartridge Systems

High Temperature Stainless Steel Cartridge System

This unique system is ideal when water system temperature routinely exceeds 100 °F (37.8 °C). Available in single or dual cartridge designs, 316 gauge stainless steel withstands higher temperatures and allows for pressures up to 160 psi.

Typical applications - cooling loops

Economy Filter Holder System

This filter holder is a basic holder system. It accepts 10 in (25.43 cm) double open-ended filters and can be piped into any water line.

Typical applications - particle filtration

Harvey DI+ Cartridge System

This turnkey system includes a remote dispenser, Pura-Lite (50 K Ω -cm) indicator and a high capacity, two-bed ion exchange cartridge to remove calcium, magnesium, and other ions from water. It accepts a broad range of other cartridges to remove impurities and produces pure water at a flow rate of 0.5 L/min.

Typical applications - pretreatment, deionization



Thermo Scientific Barnstead

Pressurized cartridge systems

Pressurized Cartridge Systems

	Inlet Feed Line NPTF inches (cm)	Feed Water Temp. °F (°C)	Dimensions W x D x H inches (cm)	Max Feed Water Pressure (psig)	Max Flow Rate (L/min)	Voltage (VAC)	Model No.
B-Pure Single Holder	1/2 (1.3)	40-120 (4-49)	7 x 7 x 15 (17.8 x 17.8 x 38.1)	100	4	120	D4511
B-Pure Double Holder with Digital Purity Meter	1/2 (1.3)	40-120 (4-49)	15 x 7 x 27 (38.1 x 17.8 x 68.6)	100	4	120	D4521
						240	D4522-33
B-Pure Double Holder with Pura-Lite Indicator (50 kΩ-cm)	1/2 (1.3)	40-120 (4-49)	15 x 7 x 27 (38.1 x 17.8 x 68.6)	100	4	120	D4524
B-Pure Double Holder with Pura-Lite Indicator (200 kΩ-cm)	1/2 (1.3)	40-120 (4-49)	15 x 7 x 27 (38.1 x 17.8 x 68.6)	100	4	120	D5831
B-Pure Double Holder with Pura-Lite Indicator (1 MΩ-cm)	1/2 (1.3)	40-120 (4-49)	15 x 7 x 27 (38.1 x 17.8 x 68.6)	100	4	120	D5833
B-Pure Half-size Holder	1/2 (1.3)	40-120 (4-49)	7 x 7 x 15 (17.8 x 17.8 x 38.1)	100	2	120	D5839
High Temperature Stainless Steel Single Holder	1/4 (.64)	40-180 (4-82)	6 x 7 x 22 (15.2 x 17.8 x 55.9)	100	3	n/a	D0807
High Temperature Stainless Steel Double Holder	1/4 (.64)	40-180 (4-82)	13 x 7 x 22 (33 x 17.8 x 55.9)	100	3	n/a	D0855
Economy Filter Holder Use with filters only	3/4 (1.9)	40-120 (4-49)	5.25 x 5.25 x 12 (13.3 x 13.3 x 30.5)	100	2	n/a	15840
Harvey DI+ Ships with Cartridge (D0803), Pura-Lite Indicator (E3450) and Dispenser	1/4 (.64)	40-120 (4-49)	7 x 13 x 21 (17.8 x 33 x 53.3)	100	4	120	AY1273X4

Optional Accessories for Pressurized Cartridge Systems

Description	Model No.	
	120V	240V
Digital Purity Meter Compatible with B-Pure Double Holders exclusively	Resistivity range of 0.1 - 18.2 MΩ-cm	
Pura-Lite Indicator Compatible with all systems above, except with Economy Filter	Resistivity range of 50 KΩ-cm	
	Resistivity range of 200 KΩ-cm	
	Resistivity range is up to 1 MΩ-cm	
Flow Meters Compatible with all systems above, installed in the feed line	1 to 75 LPH using 1/8 in (.32 cm) NPT inlet	
	0 to 190 LPH using 1/8 in (.32 cm) NPT inlet	
0.2µm Final Filter For use with B-Pure Double Holders exclusively	D3750	D3750
Dual Pressure Gauge Kit 0 to 160 psi pressure gauge exclusive to the B-Pure Half-size System. Pressure drop greater than 10 psi indicates need for filter change	D2780	D2780
Single Pressure Gauge Kit 0 to 160 psi pressure gauge exclusive to the B-Pure Half-size System. Pressure drop greater than 10 psi indicates need for filter change	D0780	D0780

Cartridges and Filters for Pressurized Systems

Cartridge	Resin Type(s)	Capacity as CaCO ₃	Model No.
Cartridges Compatible with B-Pure Double/Single Holder and Harvey DI+ Cartridge Systems			
Anion Removal	Strong Anion Exchange	1680 gr	D0760
Cation Removal	Strong Cation Exchange	3170 gr	D0815
High Capacity, Two-Bed Ion Exchange	Two-Bed Ion Exchange	1760 gr	D0803
Macroreticular and Organic Removal	Macroreticular, Activated Carbon	2000 gal	D0836
Organic Removal	Activated Carbon	5000 gal	D0813
Oxygen Removal	Oxygen Removal	30 g	D0811
Ultrapure	Mixed-Bed	875 gr	D0809
Ultrapure and Organic Removal	Mixed-Bed, Activated Carbon	785 gr/ 1000 gal	D0832
Ultrapure and Oxygen Removal	Ultrapure, Oxygen Removal	455 gr/ 30 g	D8809
Ultrapure, Oxygen, and Organic Removal	Ultrapure, Oxygen Removal, Organic Removal	365 gr/2000 gal /12 g	D8811
Cartridges Compatible with B-Pure Half-size Holder Cartridge Systems			
Organic Removal	Activated Carbon	5000 gal	D50215
Mixed-Bed with Organic Removal	Mixed-Bed Activated Carbon	275 gr/2000 gal	D50217
Oxygen Removal	Porous Strong Anion Exchange	18 g	D50214
Mixed-Bed with Oxygen Removal	Oxygen Removal Mixed-Bed	275 gr/ 9g	D50216
Ultrapure	Mixed-Bed	550 gr	D50213
Cartridges Compatible with High Temperature Stainless Steel Double/Single Holder Cartridge Systems			
High Temp, Organic Removal	Activated Carbon	5000 gal	D8815
High Temp, Oxygen Removal	Oxygen Removal	30 g	D8817
High Temp, Ultrapure	Mixed-Bed	915 gr	D8816
High Temp, Ultrapure, Oxygen, and Organic Removal	Ultrapure, Oxygen Removal, Organic Removal	365 gr / 2000 gal / 12 g	D8818
High Temp, Mixed-Bed, Oxygen Removal	Ultrapure, Oxygen Removal, Organic Removal	455 gr / 30 g	D8825
Filters Compatible with B-Pure Half-size Holder Cartridge Systems			
1 µM Pre-filter			FL583X4
5 µM Pre-filter			FL583X1
10 µM Pre-filter			FL583X2
15 µM Pre-filter			FL583X3
0.2 µM Pre-filter			FL583X6
0.45 µM Pre-filter			FL583X5
Filters Compatible with Economy Filter Holder			
5 µM Pre-filter			D2729
10 µM Pre-filter			18011
15 µM Pre-filter			18018

Capacity (grains shown as CaCO₃)

gr = grain | gal = gallon | g = grams of oxygen

CALCULATING ION EXCHANGE CAPACITY

$$\frac{\text{Cartridge Grain Capacity}}{\text{Total Dissolved Solids (Grains/Gal)}} = \text{Gallons Processed}$$

CARTRIDGE LIFE CALCULATION

$$\frac{\text{Cartridge Grain Capacity}}{[\text{Feed Water TDS (ppm)}]} \times 64.7 = \text{Pure Water Output (Liters)}$$

$$\frac{\text{Cartridge Grain Capacity}}{[\text{Feed Water TDS (ppm)}]} \times 17.1 = \text{Pure Water Output (Gallons)}$$

Thermo Scientific Barnstead Non-pressurized cartridge systems

Hose Nipple Cartridge Systems

An economical and simple approach for purifying water, this system can be installed by itself or in a series using simple tubing. Cartridges for this system provide a selection of removal properties, such as organic, ion, metal, and oxygen removal within a 100% virgin polypropylene cartridge housing. This holder offers the flexibility of accepting half-size, two-third size and full-size cartridges as needed. Cartridges are also available with color change indicators so you know when it is time to replace the cartridge.

Typical applications - deionization, organic removal, chlorine removal



Shown here installed in a series for greater capacity, our most common Hose Nipple Cartridge, D8901 and D8902, is used for generating deionized water for glassware washing and rinsing



Bantam Deionizer Cartridge Systems

A good choice for deionization or pretreatment of lab water. Ideal for applications requiring up to 38 L/hr of purified water, this system includes an integrated resistivity meter that measures up to 18 M Ω - cm tolerating inlet pressure ranging from 5 -70 psi. Designed for non-pressurized applications, this system requires diverting output water to an atmospherically vented receptacle as no back pressure can be accepted.

Typical applications - deionization



Non-pressurized Cartridge Systems

	Inlet Feed Line NPTF inches (cm)	Feed Water Temp. °F (°C)	Dimensions W x D x H inches (cm)	Max Feed Water Pressure (psig)	Max Flow Rate (L/min)	Voltage (VAC)	Model No.
Hose Nipple Cartridge Holder	n/a	40-120 (4-49)	n/a	n/a	1.25	n/a	D8900
Bantam Deionizer Cartridge System Integrated purity indicator (25 kΩ - 18 MΩ-cm)	3/8 (.95 cm)	40-120 (4-49)	6 x 9 x 28.5 (15.2 x 22.2 x 72.4)	70	0.3	120 240	D0800 D0805

Cartridges for Non-pressurized Systems

Cartridge	Resin Type(s)	Color Indicator*	Capacity as CaCO ₃	Model No.
Cartridges Compatible with Hose Nipple Cartridge Systems				
1/2 Size Mixed-Bed	Ultrapure Ion Exchange	no	430 gr	D50220
1/2 Size Mixed-Bed with Oxygen Removal**	Ultrapure Mixed-Bed with Strong Anion Exchange	no	280 gr/4.4 g	D8822
2/3 Size High Capacity**	Two-Bed Ion Exchange with Strong Cation Exchange	yes	1100 gr	D8950
2/3 Size Mixed-Bed and Organic Removal**	Mixed-Bed Ion Exchange with Activated Carbon	yes	470 gr/200 gal	D8951
Full Size Cation Removal	Cation Exchange	yes	3000 gr	D8905
Full Size High Capacity	Two-Bed Ion Exchange	yes	1650 gr	D8901
Full Size Macroreticular, Cation Removal	Macroreticular Cation Exchange	no	n/a	D8908
Full Size Organic Removal	Activated Carbon	no	5000 gal	D8904
Full Size Oxygen Removal	Strong Anion Exchange	no	30 g	D8903
Full Size Pretreatment and Scale Eliminator	Mixed-Bed and Strong Cation Exchange, Activated Charcoal	yes	1250 gr/1000 gal	D8921
Full Size Ultrapure	Ultrapure Ion Exchange	no	915 gr	D8911
Full Size Ultrapure	Ultrapure Mixed-Bed Ion Exchange	yes	915 gr	D8902
Full Size Ultrapure with Organic Removal	Mixed-Bed Ion Exchange with Activated Carbon	yes	730 gr/2000 gal	D8922
Cartridges Compatible with Bantam Deionizer Cartridge Systems				
Anion Removal	Strong Anion Exchange	no	1680 gr	D0760
Cation Removal	Strong Cation Exchange	no	3170 gr	D0815
High Capacity, Two-Bed Ion Exchange	Two-Bed Ion Exchange	no	1760 gr	D0803
Ultrapure	Mixed-Bed	no	875 gr	D0809
Ultrapure and Organic Removal	Mixed-Bed, Activated Carbon	no	785 gr/ 1000 gal	D0832

*Hose Nipple Cartridge with color indicator to notify user cartridge is expired; Alcohol-containing samples cannot be used with these models

**Hose Nipple Cartridge with tapered straight nipple for 3/8in (.95 cm) ID tubing

Capacity (grains shown as CaCO₃)

gr = grain | gal = gallon | g = grams of oxygen

CALCULATING ION EXCHANGE CAPACITY

$$\frac{\text{Cartridge Grain Capacity}}{\text{Total Dissolved Solids (Grains/Gal)}} = \text{Gallons Processed}$$

CARTRIDGE LIFE CALCULATION

$$\frac{\text{Cartridge Grain Capacity}}{[\text{Feed Water TDS (ppm)}]} \times 64.7 = \text{Pure Water Output (Liters)}$$

$$\frac{\text{Cartridge Grain Capacity}}{[\text{Feed Water TDS (ppm)}]} \times 17.1 = \text{Pure Water Output (Gallons)}$$

Common Applications for Cartridge Systems

Thermo Scientific Barnstead cartridge systems are ideal for treating water used in many common laboratory settings:

Supply purified water for general lab equipment

(CO₂ incubators, environmental chambers, and water baths)

Pure water is frequently required for humidification and decontamination in common laboratory equipment. The equipment often requires water with a resistivity ranging from 50 kΩ - 1 MΩ. The use of ultrapure water may damage equipment, yet poor quality water can cause mineral deposits or contamination. Water used with this equipment typically does not need to be sterile.



Thermo Scientific Environmental Chamber



Thermo Scientific General Purpose Water Bath




Thermo Scientific Forma Direct Heat CO₂ Incubator

Cartridge Systems Suggested for General Lab Equipment

Application	Solution	System Model No.	Cartridge Model No.
Supply water for one piece of lab equipment, such as a CO ₂ incubator or a water purification system	Hose Nipple Cartridge System with High Capacity Cartridge <ul style="list-style-type: none"> • Non-pressurized system • High capacity ion exchange cartridge produces required resistivity and features color indicator to notify user when cartridge is exhausted. 	D8900	D8901
Supply water for multiple pieces of lab equipment, such as a water bath or environmental chamber	B-Pure Double Holder Cartridge System including Pura-Lite Indicator with Pretreatment, Organic, and Ultrapure Cartridges <ul style="list-style-type: none"> • Pressurized system accommodating two cartridges for increased capacity • Pura-Lite indicates when to change cartridge • Pretreatment, organic, and an ultrapure cartridge extends useful life of cartridges while removing organics and ions 	D4524	D0832 and D0809
Supply water for a couple of pieces of lab equipment and rinse dishes	Harvey DI+ System <ul style="list-style-type: none"> • Complete out of the box, Harvey DI+ includes a single pressurized cartridge system, Pura-Lite Indicator (E3450), a high capacity cartridge (D0803), and a dispenser 	AY1273X4	included

Purify the water used in cooling loops

Impurities in water can adversely impact the operational life of laser, radar, sonar, and other heat exchange systems equipped with water cooling loops. Untreated or improperly treated cooling loop water can lead to corrosion, scaling, inefficient heat transfer, and ultimately system failures.



Cartridge Systems Suggested for Typical Water Cooling Loops

Application	Solution	System Model No.	Cartridge Model No.	Accessory Model No.
Cooling loops where cartridge system will experience back pressure; normal temperatures, is not open to atmosphere	B-Pure Single Holder Cartridge System, Ultrapure/Oxygen Removal Cartridge and Pura-Lite Indicator <ul style="list-style-type: none"> • Pressurized system to accommodate back pressure • Ultrapure/oxygen removal cartridge produces required resistivity and maintains low oxygen level to prevent oxidation • Pura-Lite indicates when to change cartridge 	D4511	D8809	E3450
Cooling loops where cartridge system will NOT experience back pressure; no extreme temperatures; open to atmosphere	Hose Nipple Cartridge System with High Capacity Cartridge <ul style="list-style-type: none"> • Non-pressurized system • High capacity ion exchange cartridge produces required resistivity and features color indicator to indicate when to change cartridge. Since the loop is open to the atmosphere, there is no need for oxygen removal 	D8900	D8901	n/a
Cooling loops where temperature is typically 120.2 °F (49 °C); cartridge system may experience back pressure; is not open to atmosphere	High Temperature Cartridge System, Mixed-Bed Oxygen Removal Cartridge and Pura-Lite Indicator <ul style="list-style-type: none"> • Non-pressurized system • High capacity ion exchange cartridge produces required resistivity. Since the loop is not open to the atmosphere, there is a need for oxygen removal 	D0855	D8816 and D8825*	E3450

* See other cartridge options on page 11

Cartridge and Filter System Applications

- **Microbiology** - nutrient media preparation; glassware rinsing; environmental monitoring
- **Biochemistry** - sample preparation; spectroscopy
- **Chemistry** - sample preparation; dilutions; reagent prep
- **Environmental sciences** - use in Biochemical Oxygen Demand (BOD); Total Soluble Solids (TSS); Mixed-Liquor Volatile Suspended Solids (MLVSS)
- **Education** - fly food preparation; humidifiers

Customized Cartridges

Don't see the right cartridge for your needs?

We can custom build cartridges. Call your local sales representative to walk you through the process.

Let our experts recommend the right water purification products for you!

Request a **FREE, no obligation** Thermo Scientific H₂O Select Analysis Kit

Our comprehensive water analysis program is designed to identify the best water system and ancillary equipment based on a sample of your feed water. Upon evaluation, you will receive a recommendation along with estimated cartridge life. Request one today at www.thermoscientific.com/select.



For more information about these products or our complete water purification portfolio, please visit www.thermoscientific.com/purewater

© 2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America: USA/Canada +1 866 984 3766 (866-9-THERMO)

Europe: Austria +43 1 801 40 0, Belgium +32 53 73 42 41, France +33 2 2803 2180, Germany national toll free 08001-536 376, Germany international +49 6184 90 6000, Italy +39 02 95059 448, Netherlands +31 76 579 55 55, Nordic/Baltic/CIS countries +358 9 329 10200, Russia +7 812 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: Australia +61 39757 4300, China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200,

Japan +81 45 453 9220, New Zealand +64 9 980 6700, Other Asian countries +852 2885 4613 **Countries not listed:** +49 6184 90 6940

WPBROCARTRIDGE 1111

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific