


# General Specification Sheet

## 3 PLY FACE MASK WITH EAR LOOPS

SKU # 16594994



SIZE	COLOR	FISHER PART
Universal	Blue	16594994

ITEM	LENGTH	ALLOWANCE
A	17.5cm	+/-1cm
B	9.5cm	+/-0.5cm
C	8cm	Minimum
D	18cm	+/-1cm

**Physical** 100% virgin spunbonded polypropylene non-linting fabric and melt-blow polypropylene non-linting fabric . Elastic band is latex free and fully encapsulated by polyester thread for added comfort

**Properties** Pattern Size: 17.5(+/-1cm) x 9.5(+/-0.5cm)

Material Weight: Outer Layer: 25(+/-5) g/m<sup>2</sup> Nonwoven Polypropylene  
 Middle Layer: 25(+/-5) g/m<sup>2</sup> Melt-blow Nonwoven fabric  
 Inner Layer: 25(+/-5) g/m<sup>2</sup> Nonwoven Polypropylene

**Size** Universal

**Color** Blue

**Packaging** 50 each packed and heat sealed polybag. 40 bags per case. 2000 each per case.

**Sterilization** Non-sterile

\* Manufactured in an ISO13485 Certified facility.

© 2017 Thermo Fisher Scientific Inc. All rights reserved.  
 Trademarks used are owned as indicated at [fishersci.com/trademarks](http://fishersci.com/trademarks).

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20  
 Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478  
 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00  
 Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303  
 Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500