

# Aspiration Basic System

Economical and compact system for comfortable operation

The Fisherbrand Aspiration Basic System is a chemical resistant oil-free aspiration station for comfortable aspiration of non-flammable chemical, biological, and medical liquids in biochemistry, microbiology and cell culture technology applications. The unit is compact and portable. It includes a diaphragm pump and a large receiving bottle with hose connectors. The system can be extended with a large range of accessories such as footswitch and handvac handpipettor.

- Superior chemical resistant diaphragm pump with long service life for chemical-duty operation
- 4L polypropylene autoclavable receiving bottle with illuminated base to aid visual level checking
- Level sensor with automatic flow stop
- Rapid coupling hose connectors with automatic isolations valves
- In-line hydrophobic filter
- 5mm ID vacuum hose for connection to your own pipettor
- Lightweight, compact and portable



## TECHNICAL INFORMATION:

Pumping speed, L/min: 15  
Ultimate pressure, mbar: 150  
Tubing supplied, mm: 5  
Dimensions l/w/h, mm: 237 x 145 x 277  
Mass, kg: 6  
Electrical supply: 90-260V, 50/60Hz with Schuko and UK plug leads  
Certification: CE marked  
Warranty: 1 year

Cat. No	Description
15567479	Aspiration Basic System, 15L/min, 150mbar, 90-260V, 50/60Hz

# Accessories



Cat. No	Description
10700854	Handvac basic pipettor kit, comprising a single channel plastic adaptor for pipettor tip with ejector, 2m silicone hose and rapid action coupling hose connector



Cat. No	Description
15572620	Handvac duo kit, comprising a single channel plastic adapter for pipettor tip with ejector, 2m silicone hose, rapid coupling hose connector and Y-connector



Cat. No	Description
15351461	Handvac professional pipettor kit, comprising a single channel plastic adapter for pipettor tip with ejector, 8-channel stainless steel adapter, 8-channel plastic adapter with ejector, table top stand and 2m silicone hose



Cat. No	Description
10412163	Eight-channel stainless steel adapter



Cat. No	Description
13594340	Eight-channel plastic adapter with ejector



Cat. No	Description
15461934	4L polypropylene receiving bottle including lid with automatic flow stop, in-line hydrophobic filter, tubing with rapid coupling hose connectors



Cat. No	Description
12967914	In-line filter



Cat. No	Description
15309804	Footswitch with ON-OFF function



Cat. No	Description
10591443	Footswitch momentary